

## BM327 Microbiology Semester 2 – Lab 1 Handout

### TASK 1A: Prepare bacterial cultures

Working in a group of **four**, you will be preparing 5 bacterial cultures, which you will need in the next laboratory session for **Error! Reference source not found.** and **Error! Reference source not found.**. As a group, you will need 12 nutrient agar plates (Protocol 1.1 You will use 8 of these plates to plate out serial dilutions of 4 different bacterial cultures (Protocol 1.2), and 4 to streak a bacterial culture for single colonies (Protocol 1.3). (Each student should prepare and plate out serial dilutions of 1 bacterial culture, and prepare one streak plate.)

#### *Protocol 1.1 : Preparation of nutrient agar plates.*

1. Label an empty Petri dish on the bottom with “nutrient agar” or “NA”. Be sure to label the bottom of the plate, and to write around the circumference (not across the centre of the plate).
2. Using aseptic technique, pour the nutrient agar plate using the molten nutrient agar provided. Pour slowly and try to avoid bubbles. Fill the plate approximately halfway with the molten agar.
3. Leave the plates untouched on your bench for ~15 min or until the agar has fully solidified. Plates will cool and solidify more quickly if placed in a single layer (not stacked).

#### *Protocol 1.2: Serial dilutions and plating for single colonies*

1. Make appropriate dilutions of the original bacterial culture to bring the anticipated concentrations of microorganisms within the range required to obtain single colonies (in this case, dilute to  $10^{-4}$ ).

Make sure to work aseptically and follow good lab practice when preparing your serial dilutions:


- Change pipette tips in between each dilution
  - Pipette carefully and accurately, making sure not to transfer any excess liquid on the pipette tip from one dilution to another
  - Mix each culture/dilution thoroughly before removing a sample, to ensure that you are collecting a representative sample. Be especially careful to mix the culture before taking the first sample, as microbial cells commonly settle to the bottom of the culture tube.
  - Label tubes carefully (writing on the side of the glass universal/bottle, not on the lid); do not mix up your dilutions.
2. Using a pipette, transfer 0.1 mL of the  **$10^{-4}$  dilution** to the agar plate. Spread the liquid evenly over the surface of a nutrient agar plate, being sure to cover the entire surface. Allow the suspension to be absorbed into the surface to the plate before inverting the plate to incubate at 37 °C for 24 - 72 hours.

3. Repeat step 2, using 0.1 mL of the  $10^{-3}$  dilution and a fresh nutrient agar plate. Note that when you are working from more-dilute to more-concentrated dilutions, you can reuse the same spreader (provided that you have not contaminated it between plates, e.g. by contact with the lab bench or other non-sterile objects.)

*Competency: Preparing a streak plate*

One of the key skills that you must develop as a microbiologist is to be able to streak microbial cultures to single colonies (Figure 2). (This is also a required skill if you wish to work in many other life science disciplines as it is essential for everyday procedures such as DNA cloning.)

To achieve single colonies on a streak plate, you must essentially perform a serial dilution using a platinum inoculating loop: on each sector of the loop, you will drag successively fewer bacterial cells across the surface of the agar. The loop is sterilised between sectors using the Bunsen burner, which kills any excess bacteria, and is then touched to the previous sector briefly to pick up a smaller number of bacteria. An isolated, individual cell will grow up to form a single colony after the plate has been incubated.



Microbiology: An Evolving Science, 4/e Figure 4.13  
Copyright © 2017 W. W. Norton & Company, Inc.

Figure 1. Streaking for single colonies.

**Before** preparing your first streak plate, use the template provided in Figure 2 to sketch a plan for how you will streak the inoculating loop across your plate. Use arrows to indicate the direction you will streak the loop, and numbers to indicate the different sectors.




Figure 2. Template for planning a streak plate.

#### Protocol 1.3: *Streaking for single colonies*

1. Flame a loop, let it cool slightly and then touch it to the bacterial colony you have chosen to streak on a new plate.
2. Pick up the bottom of the new plate (leaving the lid on the bench) with your other hand, turn it over and gently streak the loop back and forth over a small area of the plate. Use light pressure so you don't gouge the agar.
3. Return the plate to its lid so as to protect it from contaminants in the air while you re-flame the loop. Allow the loop to cool.
4. Pick up the plate again and touch the cooled loop onto an area of the agar that has not previously been streaked. (This step helps to quench the heat from the loop and so prevents it from killing any bacteria it subsequently comes into contact with.)
5. Drag the loop once through the previous streaked area and streak it out across a new area of agar. As the loop drags across the previous streak it will pick up some bacteria and distribute them across the new area of agar.
6. Repeat steps 4 - 5 two or three times as necessary to cover the entire surface of the plate (See Figure 2 and your completed streak plate template).

## TASK 1B: Plasmid profiling of strains

Working in a group of four, you will isolate plasmid DNA from (???) strains of *E. coli* (Protocol 1.4), and analyse this DNA by agarose gel electrophoresis (Protocol 1.5 – the results (photograph of your gel) will be posted for you on MyPlace.

### Protocol 1.4: Isolation of plasmid DNA

1. Transfer 1.0 ml of your *E. coli* culture to a labelled microcentrifuge tube.
2. Spin for two minutes at 14,000 rpm in a microcentrifuge, making certain that the tube is correctly balanced (check with a demonstrator).
3. Discard the supernatant carefully by pipetting and add another 1.0 ml of your *E. coli* culture. Centrifuge the cells again as in step 2.
4. Discard the supernatant carefully by pipetting and resuspend the pelleted bacterial cells in 250 µl of Buffer P1.
5. Add 250 µl of buffer P2 and gently invert the tube 4 - 6 times to mix, or until the solution becomes viscous and slightly clear. Do not allow the lysis reaction to proceed for more than 5 minutes.
6. Add 350 µl of buffer N3 and invert the tube 4 - 6 times. Centrifuge for 10 minutes at 13,000 rpm in a microcentrifuge. You must ensure that your tubes are correctly balanced before using the centrifuge – check with a demonstrator.
7. Apply the supernatant from step 5 to the **QIAprep spin column** (note: not an Eppendorf/microcentrifuge tube!!!) by decanting or pipetting. Centrifuge for 60 seconds and discard the flow-through.
8. Wash the QIAprep spin column by adding 0.5 ml Buffer PB. Centrifuge for 60 seconds and discard the flow-through.
9. Wash the QIAprep spin column by adding 0.75 ml Buffer PE. Centrifuge for 60 seconds and discard the flow-through.
10. Centrifuge for 60 seconds to remove residual wash buffer.
11. Place the QIAprep column in a clean, labelled 1.5 ml microcentrifuge tube. To elute DNA, add 50 µl Buffer EB (10 mM Tris·Cl, pH 8.5) or water to the centre of the QIAprep spin column, **let stand for one minute**, and centrifuge for one minute. The plasmid DNA should be in the eluant (you may now discard the column).

**Commented [MF1]:** Moira – I seem to remember the gel bands being really very faint for this so I thought it might be helpful to add more cells – do you have any issues with this? I can't remember what volume of culture you give them normally?

*Protocol 1.5 : Gel electrophoresis to analyse plasmid profiles.*

1. [Note: this step has been performed for you.] Prepare 1% agarose gel. Add 1g of agarose to a glass flask. Add 100 ml 1X TAE buffer. Heat agarose in a microwave till agarose is completely dissolved in the buffer. Wait for agarose to cool down to about 50°C.
2. Prepare a gel tray by taping off both ends securely and placing a comb in the desired position. Pour the agarose into the gel tray and allow this to sit undisturbed until the agarose has set.
3. When the agarose has completely set, remove the tape and the comb and place the gel into an agarose tank. [Note: ensure the gel is oriented correctly, such that the DNA samples will not run off the gel.] Fill the tank with 1X TAE buffer.
4. Mix 16 µl of your plasmid DNA with 4 µl of 5X loading buffer and load this into one well of the agarose gel. Pipet slowly and carefully, avoiding bubbles. Do not overload the well.
5. Perform electrophoresis, 100 V, until the blue indicator front is about to reach the end of the gel (usually 45 minutes).